

## SPECIFICATIONS

Item #	Height (highest setting) (armrest to floor)	Height (lowest setting) (armrest to floor)	Patient Height (highest)	Patient Height (lowest)	Weight Limit
<b>10400-1/8</b>	53"	45"	5'10"	5'2"	350 lbs.
<b>10401-1/8</b>	46"	37"	5'2"	4'6"	350 lbs.
<b>10402-1/8</b>	61"	53"	6'6"	5'10"	350 lbs.
<b>10416-1</b>	40"	31"	4'6"	4'	175 lbs.
<b>10425-8</b>	50.6"	42.6"	5'10"	4'8"	250 lbs.
<b>10426-8</b>	42.13"	36.13"	5'2"	4'6"	250 lbs.
<b>10427-8</b>	60.25"	48.25"	6'9"	5'10"	250 lbs.