

## SPECIFICATIONS

Item #	Width (between arms)	Outside Width	Width Seat	Depth Seat	Height (highest)	Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
<b>11125KD-1/2</b>	17.5"	24"	13.5"	16.5"	23"	19"	18.1 lbs.	21/42 lbs.	300 lbs.
<b>11125PSKD-1</b>	17.5"	24"	13.5"	16.5"	23"	19"	20.2 lbs.	23.5 lbs.	300 lbs.