

SPECIFICATIONS

Item #	Handle Height Highest	Handle Height Lowest	Unit Width	Unit Length	Seat (depth)	Seat (width)	Seat (height)	Product Weight	Carton Shipping Weight	Weight Limit
WINNIE LITE SUPREME	38"	32"	26.5"	24"	NA	NA	NA	9 lbs.	13 lbs.	300 lbs.
WINNIE DELUXE	38"	33"	26.5"	24"	NA	NA	NA	13 lbs.	19 lbs.	300 lbs.