

PresSsion®

COMPRESSION THERAPY

Chattanooga Group is offering a new line of PresSsion® compression therapy systems, the Multi 6 and the Multi 3. Both are truly high performance gradient sequential compression systems at an economical price. The Multi 6 is the top of the line, offering adjustable treatment parameters and therapy time control. The economical Multi 3 is designed for home therapy or clinical use.



PresSsion® Multi 6 (P.N. 4333)

- Designed for easy use.
- Use with either 3 or 6 chambered garments.
- Inflates two garments simultaneously.
- Quick inflation rate.
- Adjustable pressure with Digital Pressure Meter.
- Therapy timer from 15 to 120 minutes.
- Adjustable inflation and deflation times.
- Quiet, reliable, maintenance-free cyclone pump.

PresSsion® Multi 3 (P.N. 4331)

- Designed for easy use.
- Use with either uniform (single) or 3 chambered garments.
- Inflates two garments simultaneously.
- Quick inflation rate.
- Adjustable pressure with Digital Pressure Meter.
- Simple pre-set inflate (72 sec) and deflate (18 sec) times.
- Quiet, reliable, maintenance-free cyclone pump.

Features	PresSsion Multi 6	PresSsion Multi 3
Air Output	18 Lts/min	18 Lts/min
Pressure	30-100 mmHg	30-100 mmHg
Inflate Time	5-120 sec	72 sec
Deflate Time	5-60 sec	18 sec
Timer	15-120 min	N/A



ISO 13485 CERTIFIED

4717 Adams Road
Hixson, TN 37343 USA

1-423-870-2281

1-800-592-7329 USA

1-800-361-6661 CANADA

www.chattgroup.com

PresSsion®

COMPRESSION THERAPY

Made of strong, lightweight nylon-reinforced washable polyurethane, these garments feature a flexi-walled inflation socket with a specially designed connector to assure secure hose retention and prevent air shut-off even if the hose is kinked or flattened. Garment extenders are available from 3" to 12" (8 cm to 31 cm) wide for all garments.



PresSsion Garments (Uniform)

- 43016** – Foot and Ankle – 10" (25 cm)
- 43017** – Half Leg – 20" (51 cm)
- 43018** – 3/4 Leg – 26" (66 cm)
- 43019** – Full Leg – 29" (74 cm)
- 43020** – Full Leg – 32" (81 cm)
- 43021** – Full Leg – 35" (89 cm)
- 43023** – Full Leg Wide – 35" (89 cm)
- 43013** – Half Arm – 17" (43 cm)
- 43014** – Full Arm – 29" (74 cm)

PresSsion Sequential Garments

- 43036** – Half Leg 3 Chambers 20" (51 cm)
- 43194** – 3/4 Leg 3 Chambers 26" (66 cm)*
- 43092** – Full Leg 3 Chambers 32" (81 cm)*
- 43149** – 3/4 Leg 6 Chambers 26" (66 cm)
- 43032** – Full Leg 6 Chambers 32" (81 cm)
- 43034** – Half Arm 3 Chambers 17" (43 cm)
- 43093** – Full Arm 3 Chambers 30" (76 cm)*
- 43035** – Full Arm 6 Chambers 30" (76 cm)

* For use with Multi 3 only.

Garments Extenders Single Chamber

Half Leg

- 43132** – 3" (8 cm) width for 20" (51 cm)
- 43133** – 6" (15 cm) width for 20" (51 cm)
- 43134** – 12" (30 cm) width for 20" (51 cm)

Leg

- 43135** – 3" (8 cm) width for 26.75" (66 cm)
- 43136** – 6" (15 cm) width for 26.75" (66 cm)
- 43137** – 12" (30 cm) width for 26.75" (66 cm)

Full Leg

- 43138** – 3" (8 cm) width for 29" (74 cm)
- 43139** – 6" (15 cm) width for 29" (74 cm)
- 43140** – 12" (30 cm) width for 29" (74 cm)
- 43141** – 3" (8 cm) width for 32" (81 cm)
- 43142** – 6" (15 cm) width for 32" (81 cm)
- 43143** – 12" (30 cm) width for 32" (81 cm)
- 43144** – 3" (8 cm) width for 35" (89 cm)
- 43145** – 6" (15 cm) width for 35" (89 cm)
- 43146** – 12" (30 cm) width for 35" (89 cm)