

# SPECIFICATIONS

MODEL	CPN16FBA	CPN18FBA	CPN20FBA	CPN22FBA HD
SADMERC CODE	K0822	K0822	K0822	K0824
DRIVER WHEELS	12.5" BLACK PU TIRES	12.5" BLACK PU TIRES	12.5" BLACK PU TIRES	12.5" SOLID TIRES
CASTER WHEELS	8" BLACK PU TIRES	8" BLACK PU TIRES	8" BLACK PU TIRES	8" SOLID ARTICULATING
ANTI-TIP WHEELS	1.5" SOLID	1.5" SOLID	1.5" SOLID	1.5" SOLID
	REAR MOUNTED	REAR MOUNTED	REAR MOUNTED	REAR MOUNTED
MAXIMUM SPEED	5.4 KM/H	5.4 KM/H	5.4 KM/H	5.0 MPH
BRAKES	ELECTRO MAGNETIC	ELECTRO MAGNETIC	ELECTRO MAGNETIC	ELECTRONIC MAGNETIC
	BRAKING	BRAKING	BRAKING	REGENERATIVE
	"INTELLIGENT BRAKING"	"INTELLIGENT BRAKING"	"INTELLIGENT BRAKING"	"INTELLIGENT BRAKING"
GROUND CLEARANCE	4"	4"	4"	4"
TURNING RADIUS	30"	32"	33"	35"
OVERALL SIZE LENGTH	42.5" WIDTH	42.5" WIDTH	42.5" WIDTH	42.5" WIDTH
	24" LENGTH	26" LENGTH	28" LENGTH	29" LENGTH
DRIVE TRAIN	2-MOTOR, RWD	2-MOTOR, RWD	2-MOTOR, RWD	2-MOTOR, RWD
BATTERY (INCLUDED)	12V 36A*2	12V 36A*2	12V 36A*2	12V 36A*2
RANGE	15 MILE	15 MILE	15 MILE	UP TO 10 MILES
BATTERY CHARGER	OFF-BOARD, 4 AMP	OFF-BOARD, 4 AMP	OFF-BOARD, 4 AMP	OFF-BOARD, 4 AMP
MOTOR CONTROLLER	PG (60A)	PG (60A)	PG (60A)	PG (60A)
MOTOR	24VDC, 320W, 4600RPM	24VDC, 320W, 4600RPM	24VDC, 320W, 4600RPM	24 VDC, 420W, 4200RPM
WEIGHT CAPACITY	300 LBS.	300 LBS.	300 LBS.	400 LBS.
WEIGHTS BASE	90 LBS NOBATT	90 LBS NOBATT	93 LBS NOBATT	105 LBS. NOBATT
BATTERY WEIGHT	56 LBS BATT	56 LBS BASE	56 LBS BASE	56 LBS BASE
SEAT WIDTH	16"	18"	20"	22"
SEAT DEPTH	16"	16"	16"	18"
SEAT HEIGHT	17"	17"	17"	18"
BACK HEIGHT ADJUSTABLE	SECTIONAL, 17" - 21"	SECTIONAL, 17" - 21"	SECTIONAL, 17" - 21"	SECTIONAL, 18" - 22"
ARM STYLE	ADJUSTABLE HEIGHT	ADJUSTABLE HEIGHT	ADJUSTABLE HEIGHT	ADJUSTABLE HEIGHT
FOOT PLATE DIMENSION	7.8" X 6"	7.8" X 6"	7.8" X 6"	7.8" X 6"
FOOTRESTS	SWING-AWAY DETACHABLE	SWING-AWAY DETACHABLE	SWING-AWAY DETACHABLE	SWING-AWAY DETACHABLE
ELEVATING LEG RESTS	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
WARRANTY		LIFETIME ON FRAME		
		ONE YEAR ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS		
		SIX MONTHS ON BATTERY (FROM DATE OF INSTALLATION)		