



EndurO₂

Oxygen Concentrator



Small Size, Big Performance, Easy Maintenance

Compact Design converge in the NEW Endur



User friendly design

Design keeps portability in mind

- The original **SMALL** 5-liter concentrator
- **LIGHTWEIGHT** weighing only 41 lbs.
- Footprint is less than 1 ¼ sq ft
- Wide-set wheels for stability
- Integrated handles on top and on each side



Easy internal access

Maintenance is practically effortless

- Five-year limited warranty
- Easy access to all three filters without tools
- One way check valve prevents backflow
- Easy access to all internal components
- Nationwide warranty/service center system
- Hotline: 1-866-225-9804
 - maintenance/repair
 - technical questions
- Product training available
 - basic operations and maintenance DVD
 - indepth routine service and maintenance

Solid Engineering O₂ Oxygen Concentrator

Convenience is built right in

- Easy view instrument panel (65° angle)
 - power switch
 - reset button (circuit breaker)
 - flow settings
 - hour meter
 - indicator lights



Ergonomic control panel

Performance rivals larger models

- Outlet pressure 8.5 psi
- Maximum flow rate of 5 LPM
- Concentration is 93 +/- 3% O₂ (1-5 LPM)
- Sound level of 50 dBA

Item #	Description
HCS021	Oxygen concentrator without low-O ₂ alarm
HCS02	Oxygen concentrator with low-O ₂ alarm



Easy to replace filters

medline.com
1-800-MEDLINE



Specification	Model HCS021	Model HCS02
Equipped with OCI (Oxygen Concentrator Indicator)	No	Yes
Alarms	Power Failure High and Low Pressure Circuit Breaker Trip	Power Failure High and Low Pressure Low O ₂ Purity Circuit Breaker Trip
Dimensions	22" H x 15" W x 12"D	
Weight	41 lbs. (18.6 kg)	
Electrical	115 V AC, 60 Hz, 3.8 Amps, 420 W max. Double insulated cabinet	
Max. Outlet Pressure	8.5 +/- 1.0 psi (0.6 +/- 0.1 bar)	
Flow Range	1 - 5 L.P.M.	
Oxygen Concentration	93% +/- 3%	
Filter	Cabinet Compressor Intake Bacteria	
Sound Level	50 dBA +/- 2 dBA	
Operating System	Time Cycle / Pressure Swing	
Warranty	Five-year limited warranty	
Operating Environment	Ambient Temperature 50°F to 95°F (10°C to 35°C) Atmospheric Pressure: Sea Level to 7000 ft.	
Storage/Transportation	Temperature: -40°F to 158°F (-40°C to 70°C) Humidity: 10% to 100% (including condensation) Atmospheric Pressure: Sea Level to 7000 ft.	



LIMITED WARRANTY

MEDLINE Industries, Inc. warrants to original purchaser only the EndurO2 Oxygen Concentrator and its component parts to be free from defects in materials and workmanship from date of purchase for five (5) years. If within the warranty period the product fails to perform in accordance with product specifications, Medline will repair or replace without charge at its option the defective material or part. This warranty applies to labor for repairs performed by MEDLINE Authorized Warranty Centers (AWC) and does not cover cost of labor performed by the purchaser or user, a loaner unit, or rental costs while unit is being repaired. This warranty does not cover device failure due to accident, misuse, abuse, neglect, alteration, improper service, repair by unauthorized personnel or other defects not related to materials or workmanship. The warranty does not extend to non-durable components which are subject to normal wear and replacement. Medline will pay customary freight charges from the AWC to the purchaser location only. To make a warranty claim, the purchaser must return the concentrator or the component part, freight prepaid, to a Medline AWC. It is the responsibility of the purchaser to return the entire unit or remove the defective component part(s) and pack the unit or component part(s) to avoid shipping damage. Defective component part(s) must be returned for warranty inspection using the serial number as identification within thirty (30) days of return authorization date. Warranty registration cards must be returned within thirty (30) days of the purchase for this warranty to be effective.

For warranty service, contact an Authorized Service Center at 866-225-9804.

medline.com
1-800-MEDLINE (633-5463)

©2007 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. One Medline Place, Mundelein, IL 60060

Medline Canada 1-800-396-6996
 www.medline.ca canada@medline.com
 Not all products are currently available for sale in Canada.
 Please contact your Account Representative for additional information.

MKT207339/LIT506/10M/CTY6