



Pride® Lift Chairs
The Ultimate in Comfort & Style!

Comfort and style complement one another perfectly in Pride's Heritage Collection Lift Chairs. With an array of plush fabrics to fit any décor and a superior design that provides exceptional strength, Heritage Collection chairs deliver style and durability that will last for years to come. Plus, Pride's Heritage Collection offers numerous standard features, superb comfort, and quiet and smooth performance, all at a great value.

HERITAGE COLLECTION

HERITAGE COLLECTION MODELS



GL-58	GL-358PW	GL-358S	GL-358M	GL-358L	GL-358XL	GL-358XXL
2-Position Medium to Large 5'2" - 5'10"	3-Position Petite Wide Chaise Lounger 5'3" - Below	3-Position Petite Chaise Lounger 5'3" - Below	3-Position Medium to Large Chaise Lounger 5'4" - 6'0"	3-Position Large to Extra-Large Chaise Lounger 5'9" - 6'2"	3-Position Extra-Large Chaise Lounger 5'6" - 6'1"	2-Position Extra-Large Chaise Lounger 5'6" - 6'1"
					500 LBS. WEIGHT CAPACITY	600 LBS. WEIGHT CAPACITY

SPECIFICATIONS


	GL-58	GL-358PW	GL-358S	GL-358M	GL-358L	GL-358XL	GL-358XXL
Seat To Floor:	20"	19"	19"	20"	21"	20"	20"
Seat Depth:	20"	19"	18"	20"	22"	21"	22"
Seat Width:	20"	22"	18"	20"	22"	26"	30"
Top of Back to Seat:	24.5"	23.5"	23.5"	24.5"	25.5"	29"	29"
Distance From Wall:	11.5"	17"	17"	17"	17.5"	18"	18"
Overall Chair Width:	35"	39"	35"	35"	37"	39.5"	43"
Chair Weight:	108 lbs.	118 lbs.	118 lbs.	117 lbs.	120 lbs.	159 lbs.	168 lbs.

FABRIC SELECTION


NEW STANDARD FABRICS



OPTIONAL FABRICS



PADDED-SUEDE



MICRO-SUEDE



LEVANTE VINYL



NEW PLUSH VELVET COLLECTION



NEW WINDSOR MICRO-SUEDE COLLECTION



BENEFITS

- Engineered furniture grade laminate/hardwood frame
- Quiet & smooth lift system
- Heavy-duty lift actuator & scissor mechanism
- Standard head and arm covers
- Fire retardant (California specifications; bulletin 117)
- Standard, large dual pockets
- Removable back - no tools required
- 375 lbs. weight capacity (GL-358XL supports up to 500 lbs.; GL-358XXL supports up to 600 lbs.)
- Integrated battery backup in case of power failure¹
- Transformer features self-diagnostic electronics with 12' cord
- Lifetime warranty on steel lift frame
- 7-year pro-rated warranty

OPTIONS*

- Heat & Massage
- Head pillow
- Lumbar pillow

Additional lead-time may apply for some options.

¹. Batteries not included.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.