

SPECIFICATIONS

ITEM #	GT807 (20" SEAT)	GT808 (22" SEAT)
L X W X H	47.2" X 24" X 44.5"	47.2" X 24" X 44.5"
SADMERC CODE GROUP	K0807	K0808
FRONT WHEELS	4" X 10"	4" X 10"
REAR WHEELS	4" X 10"	4" X 10"
WEIGHT WITH BATTERIES	220 LBS.	234 LBS.
SPEED	0~6 MPH	0-6 MPH
BATTERY (INCLUDED)	U1 OR 50 AMP	U1 OR 50 AMP
GROUND CLEARANCE	3.9"	3.9"
WEIGHT CAPACITY SEAT	400 LBS.	500 LBS.
SEAT/BACK HEIGHT	20" SEAT / 17" BACK	22" SEAT / 20" BACK
TRAVEL RANGE	15-20 MILES	14-18 MILES
SAFE CLIMBING ANGLE	8° +/- 1°	8° +/- 1°
SEAT TO FOOT DECK HEIGHT	17" - 20"	17" - 20"
TURNING RADIUS	51.4"	51.4"
CHARGER	OFF BOARD 4 AMP	OFF BOARD 4 AMP
BRAKE ELECTRIC/ HAND BRAKE	STANDARD	STANDARD
CONTROLLER	CURTIS CONTROLLER 110A	CURTIS CONTROLLER 110A
WARRANTY	LIFETIME ON FRAME, ONE YEAR ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS SIX MONTHS ON BATTERY (FROM DATE OF INSTALLATION)	LIFETIME ON FRAME, ONE YEAR ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS SIX MONTHS ON BATTERY (FROM DATE OF INSTALLATION)