

SPECIFICATIONS

ITEM #	2800ECBL-RCL, 2800ECBU-RCL, 2800ECBL-RCL-20, 2800ECBU-RCL-20
L X W X H	42" X 24" X 44.5"
BASE WEIGHT	88.1 LBS.
BATTERY WEIGHT	53 LBS.
SEAT WEIGHT	51 LBS.
WEIGHT (CHAIR & BATTERY)	192.1 LBS.
FRONT WHEEL	6.0" X 2.0 SOLID
DRIVE WHEEL	10" 3.3" SOLID
ANTI TIP WHEEL	6.0" X 2.0" SOLID
BATTERY	36 AMP X 2 U1
CHARGER	3A OFF-BOARD
MAX SPEED	4.0 MPH
CRUISING RANGE	12.4 ~ 18.6 MILE
WEIGHT CAPACITY	300 LBS.
MAXIMUM CLIMBING ANGLE	8°
CONTROLLER	PG VR2 50A
SEAT	18" MID BACK CAPTAIN SEAT, 20" SEAT OPTION AVAILABLE
ARMREST	ADJUSTABLE
WARRANTY	LIFETIME ON FRAME, ONE YEAR ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS, SIX MONTHS ON BATTERY