

SPECIFICATIONS

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
MG12170	15.75"	16"	14.75"	16"	21.25"	16.5"	7.5 lbs.	10.5 lbs.	300 lbs.
MG12171	15.75"	16"	14.75"	16"	21.25"	16.5"	6 lbs.	9 lbs.	300 lbs.