

# SPECIFICATIONS

---

Item #	Handle Height (highest)	Handle Height (lowest)	Unit Width	Unit Length	Seat Depth	Seat Width	Seat Height	Product Weight	Carton Shipping Weight	Weight Limit
<b>726RD-RTL</b>	37"	32"	24"	28"	14"	14"	21"	12 lbs.	16 lbs.	300 lbs.
<b>R6 SERIES</b>	39"	35"	24"	28"	14"	14"	22"	12 lbs.	16 lbs.	300 lbs.