

The patient-preferred self-treatment system that relieves pain and accelerates post-operative recovery.

put a wrap on pain



VITALWEAR

RECOVERYWRAP[®]
SYSTEM

The new standard in cold/compression therapy: portable, patient-friendly and proven effective.

The RecoveryWrap system has the best all-around design of any cold/compression device to treat pain with unprecedented ease and flexibility. With the rechargeable battery pack, it's fully portable — for use at home, work, and on the go — with a look patients feel comfortable using anywhere. It maintains a 40° - 70°F temperature safety zone and features a clear ice-water canister that's easy to monitor and refill. This system is compatible with all the same great wraps from the VitalWrap product line such as the CPM compatible knee wraps and single hand applied shoulder wraps.

SMART:

Best performance and design of any cryotherapy device on the market

PORTABLE:

Rechargeable battery pack enables therapy when and where patients need it

CONVENIENT:

User-friendly features like a clear canister and knee wrap(s) that work in CPM machines

SAFE:

Dial-controlled temperature stays in the 40° - 70°F comfort range, so it's never too cold and always comfortable



One versatile device with wide-ranging applications:

Orthopedics and obstetrics • Chiropractic • Physical therapy • Sports medicine and athletics • Emergency medicine • and many more



VITALWEAR



put a wrap on pain

Systems

order no.

VIT-00404 | RecoveryWrap System, 15VDC
(control unit, tubing set, standard RW)

Bundles

VIT-00404-1	RW System with 6", with Battery
VIT-00404-2	RW System with Medium Knee
VIT-00404-3	RW System with Medium Knee, with Battery
VIT-00404-4	RW System with Medium Shoulder
VIT-00404-5	RW System with Medium Shoulder, with Battery
VIT-00404-6	RW System with Large Knee
VIT-00404-7	RW System with Large Knee, with Battery
VIT-00404-8	RW System with Large Shoulder
VIT-00404-9	RW System with Large Shoulder, with Battery
VIT-00404-10	RW System with Wrist
VIT-00404-11	RW System with Wrist, with Battery
VIT-00404-12	RW System with Back
VIT-00404-13	RW System with Back, with Battery
VIT-00404-14	RW System with Ankle
VIT-00404-15	RW System with Ankle, with Battery
VIT-00404-16	RW System with 6"
VIT-00404-17	RW System with Hip
VIT-00404-18	RW System with Hip, with Battery

Wraps

VIT-00255	Medium VitalWrap 4.25" x 53"	for elbow and calf
VIT-00001	Standard VitalWrap 6" x 60"	general purpose wrap for hip, calf, knee, hamstring, back, chest and abdomen
VIT-00118	Large VitalWrap 8" x 72"	broad coverage wrap for leg and back
VIT-00388	Back VitalWrap	for use on back
VIT-00119	Medium Shoulder VitalWrap	for use on either shoulder
VIT-00395	Large Shoulder VitalWrap	for use on either shoulder
VIT-00394	Wrist VitalWrap	for hand and wrist
VIT-00396	Hand VitalWrap	for arthritic hands
VIT-00397	Medium Knee VitalWrap	for post-surgery and CPM applications
VIT-00392	Full Knee VitalWrap	for post-surgery and CPM applications
VIT-00398	Ankle VitalWrap	for foot and ankle
VIT-00122	Cervical VitalWrap	for neck applications
VIT-00648	Hip VitalWrap	for use on hip — left and right

Accessories

VIT-00094	Carrying Case
VIT-00034	VitalWrap Tubing Set
VIT-00496	RW Wall Adapter Power Supply
VIT-00491	RW Battery Pack

SADMERC Durable Medical Equipment

E0218 | Water circulating cold pad with pump

SYSTEM SPECIFICATIONS:

PHYSICAL

size (approximate)	8" diameter x 11.8" height (20.3 cm x 30.0 cm)
weight (dry)	5 lbs (2.27 kg)
fluid tether	8' (2.4 m)

CONTROL SYSTEM

type	manual user control
------	---------------------

THERMAL SYSTEM

range	40° - 70° F (4.4° - 21.1° C)
contrast therapy transition time	1 min typical from extremes in temperature range
operating period	2-6 hours typical continuous use

CIRCULATING SYSTEM

reservoir capacity	0.8 gal (2.9 L)
reservoir fluid	ice water
flow rate thru VW	5.5 gph (21 Lph)
maximum pressure	18 psi

ELECTRICAL SYSTEM:

voltage	15 VDC
max power	15 W
leakage current	under 300 µA

CLASSIFICATION:

classification	class II equipment
type of equipment	type BF

00735 Rev. C

VitalWear, Inc.
384 Oyster Point Blvd., Suite 16
South San Francisco, CA 94080-1968
Tel: 1-800-553-4081
www.vitalwear.com