

SPECIFICATIONS

ITEM #	S35001692 (RED), S35002692 (BLUE)
L X W	37.8" X 21"
SADMERC CODE GROUP	GROUP 1/K0800
FRAME WEIGHT	47 LBS.
BATTERY WEIGHT	21 LBS.
SEAT WEIGHT	14 LBS.
GROSS WEIGHT (WITH/WITHOUT BATTERY)	82/61 LBS
DYNAMIC STABILITY INCLINE IN DEGREES	10°
MINIMUM OBSTACLE CLIMBING ABILITY	1"
SEAT DIMENSIONS	15.5" X 13.5"
SEAT TO FOOT DECK HEIGHT	15" - 17"
DRIVE WHEEL	7.5" SOLID
CASTER WHEEL	8" SOLID
BATTERY (INCLUDED)	12AH, 12V X 2
CHARGER	2 AMP OFF BOARD
MAX. SPEED	4 MPH
MAX. WEIGHT LOAD	250 LBS
CRUISING RANGE	10 MILES
CLIMBING ANGLE (SUGGESTED)	6°
TURNING RADIUS	31"
GROUND CLEARANCE	2"
MOTOR	270W X 1
ANTI-TIP WHEEL	2"
CHARGING TIME, HOURS	APPROX 6 HR
CONTROLLER	PG S-DRIVE 45A
WARRANTY	LIFETIME ON FRAME, ONE YEAR ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS, SIX MONTHS ON BATTERY (FROM DATE OF INSTALLATION)