

SPECIFICATIONS

ITEM #	S35001697 (RED), S35002697 (BLUE)
L X W	37.8" X 21"
SADMERC CODE GROUP	GROUP 1/K0814
SEAT TO FLOOR HEIGHT	20" - 22"
TOTAL WEIGHT WITH BATTERY	132 LBS
BATTERY WEIGHT	20 LBS
BACK WEIGHT	31 LBS
FRONT WEIGHT	24 LBS
SEAT WEIGHT	27 LBS
DRIVE WHEEL	8" SOLID TIRE
CASTER WHEEL	6" SOLID TIRE
ANTI-TIP WHEEL	2"
BATTERY (INCLUDED)	12AH/12V X 2
CHARGER	2A OFF-BOARD
TURNING RADIUS	19.7"
GROUND CLEARANCE	2.1"
MAX. SPEED	4.0 MPH
CRUISING RANGE	10 MILES (UP TO 20 MILES WHEN BATTERIES ARE STACKED)
MINIMUM OBSTACLE CLIMBING ABILITY	25 MM
WEIGHT CAPACITY	250 LBS.
MAXIMUM CLIMBING ANGLE	6°
CONTROLLER	60 AMP, PG VR-2
WARRANTY	LIFETIME ON FRAME ONE YEAR ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS SIX MONTHS ON BATTERY (FROM DATE OF INSTALLATION)