

SPECIFICATIONS

ITEM #	S35005GT, S35006GT
L X W	44.4" X 22.5"
SADMERC CODE GROUP	GROUP 2/KO806
SEAT DIMENSIONS	17"(D) X 17.5"(W)
SEAT TO FOOT DECK HEIGHT	17" - 20"
GROSS WEIGHT (WITH/WITHOUT BATTERY)	125/94 LBS
REAR WEIGHT	42 LBS
FRONT WEIGHT	26 LBS
SEAT WEIGHT	26 LBS
BATTERY WEIGHT TOTAL	31 LBS
TIRE SIZE (FRONT)	10" SOLID INSERTS
TIRE SIZE (REAR)	10" SOLID INSERTS
BATTERY (INCLUDED)	20 AMP X 2
CHARGER	3 AMP
MAX. SPEED	5 MPH
MAX WEIGHT LOAD	300 LBS
CRUISING RANGE	15 MILES (W/20 AMP BATTERY) 20 MILES (W/ U1 BATTERY)
SAFE CLIMBING ANGLE (SUGGESTED)	8°
TURNING RADIUS	40.2"
GROUND CLEARANCE	2.5"
CONTROLLER	DRIVE CTE 70 AMP
WARRANTY	LIFETIME ON FRAME ONE YEAR ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS SIX MONTHS ON BATTERY (FROM DATE OF INSTALLATION)