

SPECIFICATIONS

ITEM #	S35010
L X W	37.8" X 20"
SADMERC CODE GROUP	GROUP 1/K0800
REAR WEIGHT	24 LBS
FRONT WEIGHT	34 LBS
SEAT WEIGHT	26 LBS
BATTERY WEIGHT	20 LBS
GROSS WEIGHT	104 LBS
DYNAMIC STABILITY INCLINE IN DEGREES	10°
MINIMUM OBSTACLE CLIMBING ABILITY	1"
SEAT DIMENSIONS	17"(D) X 17.5"(W)
SEAT TO FOOT DECK HEIGHT	17-19
DRIVE WHEEL	8" SOLID
CASTER WHEEL	8" SOLID
BATTERY (INCLUDED)	12AH, 12V X 2
CHARGER	2 AMP OFF BOARD
MAX. SPEED	4 MPH
MAX. WEIGHT LOAD	250 LBS
CRUISING RANGE	UP TO 15 MILES (UP TO 30 MILES WITH STACKABLE BATTERY)
SAFE CLIMBING ANGLE	8°
TURNING RADIUS	32.3"
GROUND CLEARANCE	2"
MOTOR	270W X 1
ANTI-TIP WHEEL	2"
CHARGING TIME, HOURS	6 HR
CONTROLLER DRIVE	70A
WARRANTY	LIFETIME ON FRAME ONE YEAR ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS SIX MONTHS ON BATTERY (FROM DATE OF INSTALLATION)