

## SPECIFICATIONS

| ITEM #                               | S35015  |
|--------------------------------------|---|
| L X W                                | 37.8" X 20"   |
| SADMERC CODE GROUP                   | GROUP 1/K0800   |
| REAR WEIGHT                          | 24 LBS  |
| FRONT WEIGHT                         | 41 LBS  |
| SEAT WEIGHT                          | 26 LBS  |
| BATTERY WEIGHT                       | 20 LBS  |
| GROSS WEIGHT                         | 111 LBS   |
| DYNAMIC STABILITY INCLINE IN DEGREES | 10°   |
| MINIMUM OBSTACLE CLIMBING ABILITY    | 1"  |
| SEAT DIMENSIONS                      | 17"(D) X 17.5"(W)   |
| SEAT TO FOOT DECK HEIGHT             | 17" - 19"   |
| DRIVE WHEEL                          | 8" SOLID  |
| CASTER WHEEL                         | 8" SOLID  |
| BATTERY (INCLUDED)                   | 12AH, 12V X 2   |
| CHARGER                              | 2 AMP OFF BOARD   |
| MAX. SPEED                           | 4 MPH   |
| MAX. WEIGHT LOAD                     | 250 LBS   |
| CRUISING RANGE                       | UP TO 15 MILES (30 MILES WITH STACKABLE BATTERIES)  |
| SAFE CLIMBING ANGLE                  | 8°  |
| TURNING RADIUS                       | 44"   |
| GROUND CLEARANCE                     | 2"  |
| MOTOR                                | 270W X 1  |
| ANTI-TIP WHEEL                       | 2"  |
| CHARGING TIME, HOURS                 | APPROX 6 HR   |
| CONTROLLER                           | DRIVE 70A   |
| WARRANTY                             | LIFETIME ON FRAME<br>ONE YEAR ON ELECTRONIC CONTROLLER AND<br>DRIVE TRAIN COMPONENTS<br>SIX MONTHS ON BATTERY (FROM DATE OF INSTALLATION) |