

# SPECIFICATIONS

|                                      |  |
|--------------------------------------|--|
| ITEM #                               | S37600 (BLUE), S37601 (RED)  |
| L X W X H                            | 37.4" X 19.8" X 32.5"  |
| SADMERC CODE GROUP                   | GROUP 1/K0800  |
| REAR WEIGHT                          | 24 LBS   |
| FRONT WEIGHT                         | 26 LBS   |
| SEAT WEIGHT                          | 22 LBS   |
| BATTERY WEIGHT                       | 20 X 2 LBS   |
| GROSS WEIGHT                         | 113 LBS  |
| DYNAMIC STABILITY INCLINE IN DEGREES | 10°  |
| MINIMUM OBSTACLE CLIMBING ABILITY    | 1"   |
| SEAT DIMENSIONS                      | 16.5" (D) X 15.5." (W)   |
| DRIVE WHEEL                          | 7" SOLID   |
| CASTER WHEEL                         | 7" SOLID   |
| BATTERY (INCLUDED)                   | 12AH, 12V X 2  |
| CHARGER                              | 3 AMP OFF BOARD  |
| MAX. SPEED                           | 4 MPH  |
| MAX. WEIGHT LOAD                     | 250 LBS  |
| CRUISING RANGE                       | 10 MILES   |
| SAFE CLIMBING                        | ANGLE 8°   |
| TURNING RADIUS                       | 22.5"  |
| GROUND CLEARANCE                     | 2.5"   |
| MOTOR                                | 180W X 1   |
| ANTI-TIP WHEEL                       | 2"   |
| CHARGING TIME, HOURS                 | 6 HR   |
| WARRANTY                             | LIFETIME ON FRAME, ONE YEAR ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS<br>SIX MONTHS ON BATTERY (FROM DATE OF INSTALLATION) |