

## SPECIFICATIONS

<b>ITEM #</b>	<b>S45001GT, S45002GT</b>
L X W	49" X 22.5"
SADMERC CODE GROUP	GROUP 2/K0806
SEAT DIMENSIONS	17"(D) X 17.5"(W)
SEAT TO FOOT DECK HEIGHT	17" - 20"
GROSS WEIGHT (WITH/WITHOUT BATTERY)	136/106 LBS
REAR WEIGHT	42 LBS
FRONT WEIGHT	38 LBS
SEAT WEIGHT	26 LBS
BATTERY WEIGHT TOTAL	31 LBS
TIRE SIZE (FRONT)	10" SOLID INSERTS
TIRE SIZE (REAR)	10" SOLID INSERTS
BATTERY (INCLUDED)	20 AMP X 2
MAX. SPEED	5 MPH
MAX. WEIGHT LOAD	300 LBS
MINIMUM OBSTACLE CLIMBING ABILITY	2"
CRUISING RANGE	15 MILES (20 MILES W/OPTIONAL U1)
SAFE CLIMBING ANGLE (SUGGESTED)	8°
TURNING RADIUS	64"
GROUND CLEARANCE	3"
CONTROLLER	DRIVE 70/AMP CTE
WARRANTY	LIFETIME ON FRAME ONE YEAR ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS SIX MONTHS ON BATTERY (FROM DATE OF INSTALLATION)