

SPECIFICATIONS

ITEM #	S45200, S45300
L X W	51" X 26.4"
SADMERC CODE GROUP	K0899 (MISCELLANEOUS)
SEAT DIMENSIONS	17"(D) X 17.6"(W)
SEAT TO FOOT DECK HEIGHT	17" - 20"
GROSS WEIGHT (WITH/WITHOUT BATTERY)	201/151 LBS
REAR WEIGHT	64 LBS
FRONT WEIGHT	61 LBS
SEAT WEIGHT	26 LBS
BATTERY WEIGHT TOTAL	50 LBS
TIRE SIZE (FRONT)	13" SOLID
TIRE SIZE (REAR)	13" SOLID
BATTERY (INCLUDED)	U1 X 2
CHARGER	3 AMP
MAX. SPEED	6 MPH
MAX. WEIGHT LOAD	350 LBS
CRUISING RANGE	20 MILES
SAFE CLIMBING ANGLE	10°
TURNING RADIUS	70"
CONTROLLER DYNAMIC	75 AMP MULTI MODE
GROUND CLEARANCE	4.5"
WARRANTY	LIFETIME ON FRAME ONE YEAR ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS SIX MONTHS ON BATTERY (FROM DATE OF INSTALLATION)