

SPECIFICATIONS

ITEM #	SP-3C-701	SP-3C-701-22
POWER CHAIR SIZE	43" (L) X 26" (W) X 46" (H)	49.8" (L) X 25.1" (W) X 43.7" (H)
SADMERC CODE GROUP	K0825	K0852
WEIGHT, W/ BATTERY	219.5 LBS.	217.8 LBS
WEIGHT, W/O BATTERY	154.4 LBS.	170 LBS.
SEAT TO FLOOR HEIGHT	23" - 25"	23" - 25"
SEAT DIMENSIONS	19.5" (W) X 18.5" (D)	22" (W) X 21" (D)
BATTERY (INCLUDED)	12V 36AH X 2 OR 50 AMP	12V/36AH X 2 OR 50 AMP
MOTOR	450W 2	450W X 2 HD
CHARGER	DC24V 3A ON-BOARD	3A ON-BOARD
CASTER WHEELS	9" X 2 SOLID TIRE	9" X 2 SOLID TIRE
DRIVE WHEELS	13" X 5.5" FLAT FREE	13" X 5.5" FLAT FREE
DRIVE SYSTEM	REAR WHEEL DIRECT DRIVE	REAR WHEEL DIRECT DRIVE
BRAKE SYSTEM	ELECTROMAGNETIC BRAKES	ELECTROMAGNETIC BRAKES
TOP SPEED	5 MPH	5 MPH
CLIMBING ANGLE	8°	8°
CRUISING RANGE	14-18 MILES / 18-22 MILES (50 AMP)	12-15 MILES / 18-22 MILES (50 AMP)
TURNING RADIUS	42"	42"
ARMRESTS	ADJUSTABLE WIDTH, HEIGHT AND ANGLE	ADJUSTABLE WIDTH, HEIGHT AND ANGLE
FOOTPLATE	HEIGHT ADJUSTABLE ANGLE AND DEPTH	HEIGHT ADJUSTABLE ANGLE AND DEPTH
FOOTREST OR ELR'S	OPTIONAL	OPTIONAL
MAXIMUM WEIGHT CAPACITY	400 LBS	500 LBS
CONTROLLER	PG VR2 70A	PG VR2 70A
WARRANTY	LIFETIME ON FRAME, ONE YEAR ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS SIX MONTHS ON BATTERY (FROM DATE OF INSTALLATION)	LIFETIME ON FRAME, ONE YEAR ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS SIX MONTHS ON BATTERY (FROM DATE OF INSTALLATION)