

A revolutionary way to reduce pain and swelling through the hot, cold, contrast and compression therapies you trust.

put a wrap on pain



VITALWEAR

VITALWRAP®
SYSTEM

Did you know? 52% of Diabetics also suffer from Arthritis.

The CDC* recently reported a high incidence of arthritis among adults diagnosed with diabetes. The **SADMERC** approved VitalWrap System can help reduce pain and swelling anywhere in the body -without drugs- in one breakthrough device. With our innovative family of wraps, patients can find relief through hot, cold, contrast and compression therapies with precision and comfort. The result: chronic arthritis sufferers can get up and moving again.

SIMPLE:

Portable, easy-to-use design for patient self-treatment.

PRECISE:


Control temperature with one dial, 40-105°F.

FLEXIBLE:

Hot, cold, contrast and compression therapies.

COMPLETE:

Family of therapeutic wraps treats multiple body areas.



More gain, less pain.

52% of all Diabetics also suffer from Arthritis.

Arthritis makes it difficult for the diabetic patient to manage their condition through exercise. Now is the time to help more than 20 million people relieve joint pain and restore mobility with the SADMERC approved VitalWrap System. Register online for the full CDC report and learn how VitalWear can help diversify your product line.

Wide-ranging applications, one versatile solution.

VITALWEAR®
H0025-0000-1001 17 4 02249



* To read the CDC findings on the interaction of diabetes and arthritis, register online at:
www.vitalwear.com/cdcreport

Wide-ranging applications — One versatile solution

Arthritis • Back pain • Sprains, strains and inflammation • Post-surgical recovery
Pain management • Functional mobility • Many other conditions



VITALWEAR



put a wrap on pain

Systems

order no.

00002 VitalWrap System 120 V_{ac}, 60 Hz
(control unit, tubing set, standard VW)

00013 VitalWrap System 120 V_{ac}, 60 Hz
(control unit, tubing set, without standard VW)

Wraps

00255	Medium VitalWrap (4.25" x 52")	for elbow and calf
00001	Standard VitalWrap (6" x 60")	general purpose wrap for calf, knee, hamstring, hip, back, chest and abdomen
00118	Large VitalWrap (8" x 72")	broad coverage wrap for leg and back
00388	Back VitalWrap	for use on back
00119	Medium Shoulder VitalWrap	for use on either shoulder
00395	Large Shoulder VitalWrap	for use on either shoulder
00394	Wrist VitalWrap	for hand and wrist
00396	Hand VitalWrap	for arthritic hands
00397	Medium Knee VitalWrap	for post-surgery and CPM applications
00392	Full Knee VitalWrap	for post-surgery and CPM applications
00398	Ankle VitalWrap	for foot and ankle
00122	Cervical VitalWrap	for neck applications
00648	Hip VitalWrap	for use on hip
00967	Hip VitalWrap II	for use on hip
01061	Cap VitalWrap	for use on head

VitalPad[®] made exclusively for XBACK[®] Solution LSO Brace

01311	VitalPad for XBACK Solution, S	for XBACK Solution Model T102
01315	VitalPad for XBACK Solution, M	for XBACK Solution Model T102
01319	VitalPad for XBACK Solution, L	for XBACK Solution Model T102
01323	VitalPad for XBACK Solution, XL	for XBACK Solution Model T102
01327	VitalPad for XBACK Solution, 2XL	for XBACK Solution Model T102
01331	VitalPad for XBACK Solution, 3XL	for XBACK Solution Model T102
01335	VitalPad for XBACK Solution, 4XL	for XBACK Solution Model T102

XBACK[®] Bracing Services, Inc. is registered trademark.

Accessories

00094	Carrying Case	for easy portability of the VW System
00034	VitalWrap Tubing Set	
00987	VitalWrap "Y" Tubing Set	for use on 2 limbs
00262	Hospital Grade Power Cord	

SADMERC Durable Medical Equipment Codes

code

E0217 Water circulating heat pad with pump

E0218 Water circulating cold pad with pump

E0249 Pad for water circulation heat unit

SYSTEM SPECIFICATIONS:

PHYSICAL	
size (approximate)	8" diameter x 13.5" height (20.3 cm x 34.3 cm)
weight (dry)	5 lbs (2.27 kg)
fluid tether	8' (2.44 m)
CONTROL SYSTEM	
type	manual user control
THERMAL SYSTEM	
range	40° - 105° F (4.4° - 40.6° C)
contrast therapy transition time	1 min typical from extremes in temperature range
operating period	2-6 hours typical continuous use
thermal cutout	116.6° F (47° C)
CIRCULATING SYSTEM	
reservoir capacity	1.1 gal (4.2 L)
reservoir fluid	ice water
flow rate thru VW	5.5 gph (21 Lph)
maximum pressure	18 psi
ELECTRICAL SYSTEM:	
voltage	120 V _{ac} , 60 Hz
power	300 W
leakage current	under 300 μA
CLASSIFICATION:	
classification	class I equipment
type of equipment	type BF

01307 Rev. A

VitalWear, Inc.

384 Oyster Point Blvd., Suite 16
South San Francisco, CA 94080-1968
Tel: 1-800-553-4081
www.vitalwear.com