

ITEM #	S37650 (BLUE), S37651 (RED)
L X W X H	41" X 19.8" X 32.5"
SADMERC CODE GROUP	GROUP 1/KO800
REAR WEIGHT	24 LBS
FRONT WEIGHT	30 LBS
SEAT WEIGHT	22 LBS
BATTERY WEIGHT	20 X 2 LBS
GROSS WEIGHT	118 LBS
DYNAMIC STABILITY INCLINE IN DEGREES	10°
MINIMUM OBSTACLE CLIMBING ABILITY	1"
SEAT DIMENSIONS	16.5" (D) X 15.5." (W)
DRIVE WHEEL	7" SOLID
CASTER WHEEL	7" SOLID
BATTERY (INCLUDED)	12AH, 12V X 2
CHARGER	3 AMP OFF BOARD
MAX. SPEED	4 MPH
MAX. WEIGHT LOAD	250 LBS
CRUISING RANGE	10 MILES
SAFE CLIMBING	ANGLE 8°
TURNING RADIUS	37"
GROUND CLEARANCE	2.5"
MOTOR	180W X 1
ANTI-TIP WHEEL	2"
CHARGING TIME, HOURS	6 HR
WARRANTY	LIFETIME ON FRAME, ONE YEAR ON ELECTRONIC CONTROLLER AND DRIVE TRAIN COMPONENTS SIX MONTHS ON BATTERY (FROM DATE OF INSTALLATION)