


# FlexiFit™ 407

PREMIUM NASAL MASK

## Specification Sheet


---

A Nasal Delivery Mask




*Incorporating FlexiFit™ Technology with auto-contouring properties, the FlexiFit™ 407 Nasal Mask comes with a unique Glider™ Strap mechanism for maximum comfort and stability. Other unique features ensure a premium sleep experience.*

# Nasal Delivery - *FlexiFit™* 407




*The FlexiFit™407 Premium Nasal Mask offers an over-the-nose mask solution that fits better, feels lighter and is easier to use. The Stretchgear™ Headgear provides even greater comfort and a maintenance-free bias diffuser ensures quieter operation. The FlexiFit™407 Nasal Mask comes as an Out-of-Box Solution fully assembled and ready for use.*

## UNIQUE MASK FEATURES FOR GREATER FREEDOM, COMFORT AND COMPLIANCE




### FlexiFit™ Technology Auto-contouring fit

*With Foam Cushion*  
FlexiFit™ Technology automatically adjusts to a wide range of nasal and facial contours with simple one-step ease of fitting and no complicated adjustments.



### Glider™ Strap Anti-leak system

The sliding Glider™ Strap provides freedom of movement while maintaining a seal. Mask force remains equal when the head rotates from side to side, minimizing the occurrence of leaks.




### Stretchgear™ Headgear Deluxe mask suspension system

This is a special elastic headgear segment that maintains an equal mask force when the head tilts up or down.




### Premium Frosted Silicone Seal Optimizing seal technology

This seal technology has a soft, smooth feel to slide onto the face, and its frosted texture optimizes the mask seal.




### Efficient Bias Diffuser Quiet usage

*Maintenance Free*  
This advanced air diffuser system prevents air from blowing onto the sleeping partner and decreases noise that could potentially disrupt sleep.



### Glider™ Clip One-step removal

Simple one-step removal of headgear. No ongoing headgear adjustments required.



### Hidden Strap Release One-step removal

Easy one-step removal of headgear. Enables reassembly without having to make fitting adjustments.

## CONTENTS

- Complete mask with headgear assembled and ready for use
- Single one size Silicone Seal designed to fit most patients
- User friendly Fitting Guide

REORDER CODE: HC407

FOR MORE INFORMATION, PLEASE CONTACT YOUR LOCAL FISHER & PAYKEL HEALTHCARE REPRESENTATIVE:

**Manufacturer** Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348, Panmure, Auckland 1741, New Zealand  
Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.com Web: www.fphcare.com

**Australia** Tel: +61 3 9879 5022 Fax: +61 3 9879 5232 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **China** Tel: +86 20 3205 3486 Fax: +86 20 3205 2132 **France** Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 **Germany** Tel: +49 7182 93777 0 Fax: +49 7182 93777 99 **India** Tel: +91 80 4123 6041 Fax: +91 80 4123 6044 **Irish Republic** Tel: 1800 409 011 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 3661 7205 Fax: +81 3 3661 7206 **Middle East** Tel: +9 71 4 360 0849 Fax: +9 71 4 368 8110 **Northern Ireland** Tel: 0800 132 189 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **UK (EU Authorised Representative)** Fisher & Paykel Healthcare UK Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: +44 1628 626 136 Fax: +44 1628 626 146 **USA** Tel: 1 800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

© Fisher & Paykel Healthcare Ltd